



FOCUS ON PRECISION

FOCUS ON INNOVATION

FOCUS ON QUALITY



FOCUS ON

**HIGH SPEED TRANSFORMER**

## HIGH SPEED TRANSFORMER

THE D.O.R.C. VITRECTOMY TRANSFORMER UPGRADES YOUR STANDARD VITRECTOMY MACHINE (EG; MILLENIUM®, MEGATRON®, PENTASYS®) INTO A HIGH SPEED SYSTEM INCREASING UP TO 3 TIMES THE CUT RATE. HIGH SPEED CUTTING WITH D.O.R.C.'S LATEST 3RD. GENERATION VITRECTOMES WILL MINIMIZE TRACTION, AND IN COMBINATION WITH AN EFFICIENT ASPIRATION FLOW, OFFER EXCELLENT CUTTING PERFORMANCE. USING YOUR ORIGINAL FOOT PEDAL, THE D.O.R.C. HIGH SPEED TRANSFORMER MAINTAINS LINEAR CONTROL OVER THE MOTHER SYSTEM AND PROMOTES A SAFE VITRECTOMY IN 20, 23 AND 25 GAUGE SURGERY.



1265.HS

### 3<sup>rd</sup> GENERATION 23 GAUGE VITRECTOME

THE D.O.R.C. 23 GAUGE HIGH SPEED CUTTER PRODUCES EXCELLENT ASPIRATION FLOW IN COMBINATION WITH TRACTIONLESS CUTTING PERFORMANCE.

- CLOSEST ASPIRATION PORT TO THE TIP
- INCREASED ASPIRATION FLOW
- INCREASED SHAFT STIFFNESS
- 2500 CUTS/MIN.

2267.NMD06



3RD GENERATION HIGH FLOW AND HIGH SPEED VITRECTOME, MAX 2500 CUTS/MIN., 23G/0.6 mm (BOX/1, STERILE)

## ACCESSORIES

### VITRECTOMES – 20 GAUGE

1267.NMD25



PNEUMATIC HIGH SPEED CUTTER, MAX 2500 CUTS/MIN.,  
20G/0.9 mm (BOX/1, STERILE)

### VITRECTOMES – 25 GAUGE

1267.NMD05



PNEUMATIC HIGH SPEED CUTTER, MAX 1500 CUTS/MIN.,  
25G /0.5 mm (BOX/1, STERILE) REQUIRES 1366 ADAPTOR

### CLAES ILLUMINATED VITRECTOME

2267.IL



DISPOSABLE ILLUMINATED VITRECTOME. MAX 2500 CUTS/MIN. 20 GAUGE/0.9 mm (BOX/1, STERILE)

### POSTERIOR VITRECTOMY PACK

6279.MIL



POSTERIOR VITRECTOMY  
PACKS (BOX/5, STERILE)

CONSISTING OF:

- CARTRIDGE
- TUBING SET REF:  
1250.ABD (1X)
- 4 WAY STOPCOCK (2X)
- GPC INFUSIONSET
- END CAP (2X)
- SCREEN DRAPE
- MAYO COVER

### ACCESSORIES FOR HIGH SPEED TRANSFORMER

1265.HSA1



HIGH SPEED TRANSFORMER ADAPTOR,  
STYLE 1 "SNAP-TITE®"

1265.HSA2



HIGH SPEED TRANSFORMER ADAPTOR,  
STYLE 2 "SCHRÄDER-BELLOWS®"

1265-HSA3

HIGH SPEED TRANSFORMER ADAPTOR,  
STYLE 3, FOR GEUDER® SYSTEMS

1265-HSA4

HIGH SPEED TRANSFORMER ADAPTOR,  
STYLE 4, FOR RUCK® SYSTEMS

# FOCUS ON KNOWLEDGE

# FOCUS ON DESIGN



For additional information please contact: **D.O.R.C. International B.V.**  
Scheijdelweg 2  
3214 VN Zuidland  
The Netherlands  
Phone: +31 181 45 80 80  
Fax: +31 181 45 80 90  
e-mail: sales@dorc.nl

In US / Canada please contact: **Dutch Ophthalmic USA**  
10 Continental Drive  
Bldg 1, Exeter, NH 03833, U.S.A.  
Phone: +1 800-75-DUTCH or +1 603-778-6929  
Fax: +1 603-778-0911  
e-mail: sales@dutchophthalmicusa.com

In Latin America please contact: **D.O.R.C. Latin America**  
Scheijdelweg 2  
3214 VN Zuidland  
The Netherlands  
Phone: +1 305-663-5662  
Fax: +1 305-663-7343  
e-mail: la-sales@dorc.nl


In Scandinavia please contact: **D.O.R.C. Scandinavia AB**  
Bettorpsgatan 10  
SE-701 11 Örebro  
Sweden  
Phone: +46 19121514  
Fax: +46 19188204  
e-mail: order@dorc.se

In France please contact: **D.O.R.C. France S.A.R.L.**  
30, Boulevard Pasteur  
75015 Paris  
France  
Phone: +33 1 53 58 14 73  
Fax: +33 1 53 58 14 75  
e-mail: ventes@dorc.fr

In Germany please contact: **D.O.R.C. Deutschland GmbH**  
Charlottenstraße 68  
10117 Berlin  
Deutschland  
Phone: +49 30 201 88364  
Fax: +49 30 201 88365  
e-mail: bestellung@dorc.de



AUTHORIZED D.O.R.C. DEALER



 FOCUS MORE AT [WWW.DORC.NL](http://WWW.DORC.NL)

