



FOCUS ON PRECISION

FOCUS ON INNOVATION

FOCUS ON QUALITY



FOCUS ON  
**MICROSURGICAL SYSTEMS**



## INTRAOCULAR ELECTROMAGNET & ACCESSORIES

### INTRAOCULAR ELECTROMAGNET

1690



- CONTROLLED REMOVAL OF FERROUS INTRAOCULAR FOREIGN BODIES (IOFB)
- TRACTIVE FORCE ADJUSTABLE TO 30 GRAMS
- MAGNETIC FIELD CAN BE INITIATED OR TERMINATED VIA FOOTSWITCH CONTROL
- FOR USE THROUGH PARS PLANA, LIMBAL OR OPEN SKY SURGICAL TECHNIQUE
- SYSTEM INCLUDES ONE 1690.A HANDPIECE

### ACCESSORIES FOR INTRAOCULAR ELECTROMAGNET

1690.A.0.9

REPLACEMENT INTRAOCULAR HANDPIECE  
(20 GAUGE/0.9 mm DIAMETER)



## MICRO DIATHERMY SYSTEM & ACCESSORIES

### MICRO DIATHERMY SYSTEM

1100



- SOLID STATE, HIGH FREQUENCY DIATHERMY
- CONCENTRATES COAGULATION EFFECT AT TARGET TISSUE, MINIMIZING THERMAL DAMAGE TO ADJACENT TISSUE
- ELECTRONIC FEEDBACK CONTROL
- MONOPOLAR AND BIPOLAR ACCESSORIES ACCOMMODATE INTRAOCULAR AND EXTRAOCULAR APPLICATIONS

# MICRO DIATHERMY SYSTEM & ACCESSORIES

## MONOPOLAR DIATHERMY ACCESSORIES

**1101**  
ELECTRODE HANDLE WITH CABLE

**1102**  
ELECTRODE, WHITE, 2 mm FINE TIP; FOR PERFORATION AND RELEASE OF SUBRETINAL FLUID

**1103**  
ELECTRODE, BLACK, FOR SURFACE OR PARTIALLY PENETRATING DIATHERMY

**1104**  
ELECTRODE, GREEN, FOR SURFACE DIATHERMY

## BIPOLAR INTRAOCULAR ACCESSORIES

**1281.D206**



BACKFLUSH INSTRUMENT WITH DIATHERMY, INCL. ACTIVE ASPIRATION, 23G/0.6 mm  
ADAPTORS AVAILABLE FOR MOST VITRECTOMY SYSTEMS

**1281.D206**



**1120**



REUSABLE INTRAOCULAR COAXIAL ELECTRODE, 20 GAUGE (0.9 mm), WITH HANDLE FOR CONTROLLING RETINAL BLEEDERS

**1120.A**

REUSABLE SILICONE CABLE FOR COAXIAL DIATHERMY ELECTRODES (FOR 1120, 1121 AND 1122)

## BIPOLAR EXTRAOCULAR ACCESSORIES

**1115**



**1115**  
BIPOLAR COAGULATION FORCEPS, ANGLED

**1116**  
BIPOLAR COAGULATION FORCEPS, STRAIGHT

**1115.A**  
DIATHERMY CABLE FOR 1115/1116

**1121/1122**



**1121/1122**  
EXO-DIATHERMY HANDLE WITH CURVED OR STRAIGHT TIP, 1.3 mm

**1120.A**  
DIATHERMY CABLE FOR 1120/1121/1122

# FOCUS ON KNOWLEDGE

# FOCUS ON DESIGN



For additional information please contact: **D.O.R.C. International B.V.**  
Scheijdelweg 2, 3214 VN Zuidland  
The Netherlands  
Phone: +31 181 45 80 80  
Fax: +31 181 45 80 90  
e-mail: sales@dorc.nl

In US / Canada please contact: **Dutch Ophthalmic USA**  
10 Continental Drive  
Bldg 1, Exeter, NH 03833, U.S.A.  
Phone: +1 800-75-DUTCH or +1 603-778-6929  
Fax: +1 603-778-0911  
e-mail: sales@dutchophthalmicusa.com

In Latin America please contact: **D.O.R.C. Latin America**  
Scheijdelweg 2, 3214 VN Zuidland  
The Netherlands  
Phone: +1 305-663-5662  
Fax: +1 305-663-7343  
e-mail: la-sales@dorc.nl

In Middle East please contact: **D.O.R.C. Middle East, Africa and India**  
PO BOX 212218, Dubai, UAE  
Phone: ++971 (0) 50 911 65 87  
Fax: ++31 (0) 181 45 80 90  
e-mail: cvds@dorc.nl


In Scandinavia please contact: **D.O.R.C. Scandinavia AB**  
Bettorpsgatan 10, SE-701 11 Örebro  
Sweden  
Phone: +46 19121514  
Fax: +46 19188204  
e-mail: order@dorc.se

In France please contact: **D.O.R.C. France S.A.R.L.**  
30, Boulevard Pasteur, 75015 Paris  
France  
Phone: +33 1 53 58 14 73  
Fax: +33 1 53 58 14 75  
e-mail: ventes@dorc.fr

In Germany please contact: **D.O.R.C. Deutschland GmbH**  
Charlottenstraße 68  
10117 Berlin  
Deutschland  
Phone: +49 30 201 88364  
Fax: +49 30 201 88365  
e-mail: bestellung@dorc.de



AUTHORIZED D.O.R.C. DEALER



 FOCUS MORE AT [WWW.DORC.NL](http://WWW.DORC.NL)

